XP Dynamis

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Komple Anti-Aging Tedavileri

Estetik;

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[May 2010]

Laser and Health Academy
Review

Laser and Health Academy Review
Vol. 2010, May; www.laserandhealth.com

The Dynamis Line Nd:YAG Lasers

Content:

I. INTRODUCING THE DYNAMIS LINE OF Nd:YAG LASERS.....	30
Key features.....	30
Key documents.....	30
Technical specifications.....	30
II. Nd:YAG THE LASER OF CHOICE FOR HEATING TARGETS IN DEEPER SKIN LAYERS.....	31
Laser tissue interaction and thermal relaxation time.....	31
Nd:YAG compared to other wavelengths.....	32
Dynamis Nd:YAG laser.....	35
III. FOTONA TECHNOLOGY.....	34
Fotona VSP (Variable Square Pulse) Technology.....	34
Fotona EFC (Energy Feedback Control) Technology.....	34
Fotona SOE (Scanner Optimized Efficacy) Technology.....	35
IV. PULSE OPERATION MODE IN PRACTICE.....	35
R31 – Gaussian beam profile handpiece.....	35
R33 – Top hat profile handpiece.....	35
Pulse mode pre-sets.....	36
Pulse mode Nd:YAG clinical experience.....	37
V. SCANNER OPERATION MODE IN PRACTICE.....	38
S11 – Scanner for precision and velocity.....	38
Scanner mode pre-sets.....	40
VI. QCW OPERATION MODE IN PRACTICE.....	44
R27 – QCW handpiece.....	44
QCW clinical experience.....	45
VII. CONCLUSIONS.....	46
VIII. REFERENCES.....	47

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


- ➔ **Novel Laser Therapy in Treatment of Onychomycosis**
Jasmina Kozarev, Zdenko Vizintin
- High Fluence, High Beam Quality Q-Switched Nd:YAG Laser with Optoflex Delivery System for Treating Benign Pigmented Lesions and Tattoos**
Matjaz Lukac, Zdenko Vizintin, Marko Kazic, Tom Sult
- ➔ **Saphenous Vein Occlusion by Endovenous Laser Ablation (EVLA) with a 1064 nm VSP Nd:YAG Laser**
Andrej Sikovec
- ➔ **QCW Pulsed Nd:YAG 1064 nm Laser Lipolysis**
Matjaz Lukac, Zdenko Vizintin, Janez Zabkar, Samo Pirnat
- ➔ **Beyond Customary Paradigm: FRAC3® Nd:YAG Laser Hair Removal**
Matjaz Lukac, Martin Gorjan, Janez Zabkar, Ladislav Grad, Zdenko Vizintin
- ➔ **Parameters for the New FRAC3® Nd:YAG Laser Skin Treatment Modality**
Matjaz Lukac, Tom Sult, Janez Zabkar, Martin Gorjan, Zdenko Vizintin
- Ablation and Thermal Depths in VSP Er:YAG Laser Skin Resurfacing**
Matjaz Lukac, Tadej Perhavec, Karolj Nemes, Uros Ahcan

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
a) Long duration, low power laser pulse


b) Short duration, high power FRAC3 laser pulse

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

Spot size	3 mm	6 mm	9 mm	12 mm	15 mm	20 mm
Effective fluence relative to the effective fluence at 9 mm	62%	88%	100%	107%	112%	117%

Ref: Lukac et al, Beyond Customary Paradigm : FRAC3 ® Nd:YAG Laser Hair Removal, J of LAHA 2010, 1, 35-46.

One family, together

Hasta konforu faktörü

One family, together

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Hair removal treatment regime	Fluence	Pulse duration	Spotsize	Hair removal efficacy	Pain
Standard	55 J/cm ²	15-25 ms	9 mm	1	2
FRAC3®	55 J/cm ²	1 ms	3 mm	2.2	1


Ref: Lukac et al, Beyond Customary Paradigm : FRAC3 ® Nd:YAG Laser Hair Removal, J of LAHA 2010, 1, 35-46.

One family, together

Klinik rehberler – Laser Epilasyon

Treatment Modality:		<i>FRAC3® Hair Removal</i>		
Laser System		Fotona Dynamis		
Laser wavelength		Nd:YAG (1064 nm)		
User Interface Mode:		Pulse (Accelera)		
Utilized handpieces:		S11 scanner		
	Spot size (mm)	Pulsewidth (ms)	Fluence (J/cm²)	Frequency (Hz)
Skin group A	3	0.3- 5.0 msec	<55	15
Skin group B	3	0.3 - 5.0 msec	<35	20
Skin group C	3	0.3 - 5.0 msec	<25	25

One family, together

Nd:YAG Operasyon Modlari

- **PULSE** modu
- **SCANNER** modu
- **QCW** modu

One family, together

The Highest Performance, Best Made Laser Systems in the World

MEMORIES

Handpiece

Fluence/Energy

Frequency

>M NdPULSE PK
>R33 FRAC3
>20 J/cm2 2 mm
>1.0 Hz 0.1 ns

PERSONAL

Frac3 Gösterges

Spot size

Duration

STORE

BACK

Standby

Ready

Fotona Medical Lasers

NEW

One family, together

The Highest Performance, Best Made Laser Systems in the World

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